

# Technical Bulletin #159: Working with RNA

## Living with RNase

Most researchers are acutely aware of the risk of RNase contamination, and we do not want to belabor this point or cause undue worry. We do not routinely find microcentrifuge tubes used with RNA if they are from unopened bags or from bags in which care was taken to avoid contaminating the tubes. Yet we do consider (even those marketed as being RNase-free), the use of which results in RNA degradation. We do recommend that gloves be worn when handling any reagents which have touched refrigerator handles, door knobs, or pipettors are not RNase-free.) When performing procedures that use RNases (eg. ribonuclease protection assays), care should be taken that pipettors are not contaminated by accident. One potential source of contamination is the metal tip ejector mechanism on the side of the pipette bar when it is necessary to insert the pipettor into a larger vessel where the ejector could come into contact with the walls or contents of the vessel will eliminate

### A. Detecting RNase

While contaminating RNase can result in a failed experiment, it is often difficult and time-consuming to determine which solution or piece of equipment is responsible. No. 1964) allows researchers to identify contaminated reagents and equipment quickly, and nonisotopically. In the RNaseAlert Kit procedure, an optimized RNA substrate, containing both fluorescent and quenching moieties, is introduced as a target for any contaminating RNase. In the presence of RNase, the substrate is cleaved, releasing a fluorescence signal can be detected by eye or with a fluorometer.

### B. Getting rid of RNase

If RNase contamination of reagents or equipment is suspected to be a problem, extra precautions may be necessary. Autoclaving tips, tubes and solutions is not sufficient. Glassware can be baked at 300°C for four hours and plasticware, tubes and most solutions can be DEPC-treated (see below). However, both procedures are time-consuming and expensive and possibly carcinogenic. As an alternative, Ambion's **RNaseZap™** (Cat. No. 9780) can be used to eliminate RNase from glassware, plastic surfaces and solutions. **Zap™** has been shown to effectively inactivate 5 µg of RNase dried onto the bottom of eppendorf tubes without inhibiting subsequent enzymatic reactions performed. **Zap™** contains three ingredients known to be active against RNase. **RNase Zap™** can be poured onto or wiped over surfaces and works immediately upon contact. It is compatible with distilled water and is ready for use.

## Treating Solutions with DEPC to Remove RNase

To ensure that solutions are free of RNase contamination, they can be treated with diethylpyrocarbonate (DEPC) [WARNING: DEPC is a suspected carcinogen; handling; e.g., always wear gloves and handle under an approved fume hood]. DEPC reacts with histidine residues of proteins and will inactivate RNases. However, DEPC needs to be removed by heat treatment before the solution is used (DEPC breaks down to CO<sub>2</sub> and ethanol). Add DEPC to solutions at a concentration of 0.05% (v/v) (0.5 ml/liter); stir or shake into solution, incubate for several hours; autoclave at least 45 minutes, or until DEPC scent is gone. Please be aware that compounds containing primary amines (e.g., 2-Amino-2-hydroxymethyl-1,3-propanediol), will also react with DEPC, and thus should be added only after DEPC treatment is complete. Note: We have observed that DEPC and thoroughly autoclaved, caused a 20% inhibition of translation in a reticulocyte lysate. We find that distilled water is generally already RNase-free, and

## How to Store RNA

RNA may be stored in a number of ways. For short-term storage, RNase-free H<sub>2</sub>O (with **0.1 mM EDTA**) or **TE buffer** (10 mM Tris, 1mM EDTA) may be used. RNA can be stored for up to a year without degradation. Magnesium and other metals catalyze non-specific cleavages in RNA, and so should be chelated by the addition of EDTA if RNA is to be stored. It is important to use an EDTA solution known to be RNase-free for this purpose (older EDTA solutions may have microbial growth which could contaminate the RNA). We have suggested that RNA solubilized in formamide may be stored at -20°C without degradation for at least one year (Chomczynski, 1992).

For long term storage, RNA samples may also be stored at -20°C as ethanol precipitates. Accessing these samples on a routine basis can be a nuisance, however, and so should be pelleted and dissolved in an aqueous buffer before pipetting, if accurate quantitation is important. An alternative is to pipet directly out of an ethanol precipitate.

even suspension. We have found, however, that while this method is suitable for qualitative work, it is too imprecise for use in quantitative experiments. RNA does not precipitate well, probably because it forms aggregates; non-uniform suspension, in turn, leads to inconsistency in the amount of RNA removed when equal volumes are pipetted.

## How to Precipitate RNA

### A. Precipitating with alcohol

Precipitating RNA with alcohol (ethanol or isopropanol) requires a minimum concentration of monovalent cations (for example: 0.2 M Na<sup>+</sup>, K<sup>+</sup>; 0.5 M NH<sub>4</sub><sup>+</sup>) (Volume of alcohol added). If the cation concentration has been adjusted, the RNA may be precipitated by adding 2.5 volumes of ethanol or 1 volume of isopropanol and mixing thoroughly, followed by centrifugation. While isopropanol is somewhat less efficient at precipitating RNA, isopropanol in the presence of NH<sub>4</sub><sup>+</sup> is better than ethanol at keeping free nucleotides in the supernatant. RNA precipitation is faster and more complete at higher RNA concentrations. A general rule of thumb is that RNA concentrations of 10 µg/ml or higher will precipitate within a few hours to overnight with no difficulty, but at lower concentrations a carrier nucleic acid or glycogen should be added to facilitate precipitation and maximize recovery.

### B. Precipitating with lithium chloride

**Lithium Chloride** may also be used to precipitate RNA, and has the advantage of not precipitating carbohydrate, protein or DNA. LiCl is frequently used to remove poly(A) tails from RNA prepared by other methods. A final LiCl concentration of 2-3 M is needed to precipitate RNA (adding an equal volume of 4 M LiCl, 20 mM Tris-HCl, pH 7.5, and 100 mM NaCl will). Note that no alcohol is needed for LiCl precipitation. RNA should be allowed to precipitate at -20°C; precipitation time depends on RNA concentration. It may take several hours to overnight. After centrifugation to collect the RNA, pellets can be rinsed with 70% ethanol to remove traces of LiCl. LiCl efficiently precipitates RNA of any length. While LiCl can effectively precipitate RNA from more dilute solutions, for best results, the RNA concentration should exceed 200 µg/ml.

## Incorporation and Yield

"Percent incorporation" is calculated by comparing the amount of radioactivity incorporated into synthesized RNA with the total amount of radioactivity in the reaction mixture (see below) but can also be done by simply counting an aliquot of the transcription reaction before and after removal of unincorporated nucleotides. The comparison must be adjusted to represent equivalent aliquots. Unincorporated nucleotides may be removed by precipitation using LiCl or NH<sub>4</sub>OAc and EtOH (see below) or by reaction over an RNase-free Sephadex column (e.g., Ambion's **NucAway™ column**), or by gel purification.

The amount of radioactivity incorporated into RNA may also be determined by precipitation with trichloroacetic acid (TCA), filtration, and counting in a liquid scintillation counter. Add 98 µl of water containing 10 µg of carrier DNA or RNA. To this add 2 ml of cold 10% TCA, vortex and incubate on ice 5 minutes. Collect the precipitate by vacuum through GF/C glass fiber filters. Wash the sample tube twice with 2 ml 10% TCA and once with 2 ml of 95% ethanol, passing the washes through the filter. The filter is placed in vials with liquid scintillation cocktail and counted. Note: Both RNA and DNA may be precipitated using this method.

Since percent incorporation of a radiolabeled nucleotide is directly proportional to yield, the actual yield of a transcription reaction is equivalent to that proportion. For example, Ambion's **MAXIscript™ Kit** reactions have a theoretical 100% yield of 77 ng when the transcription reaction contains a limiting nucleotide concentration. If in a particular reaction the percent incorporation was 80%, then 0.80 X 77 ng or 62 ng of labeled RNA were synthesized.

Some ribosomal subunit size relationships within the eukaryotes are illustrated in Table 1. Both 18S and 28S rRNA contain modified nucleotides, including methylated nucleotides (e.g., 26S and 37 for 18S; 71 and 60 for 28S, respectively).

	Avg. # of bases
Organism	18S
Drosophila	1976
Rat	1874
Human	1868

Table 1. Ribosomal Subunit Sizes in Representative Eukaryotes.

## RNA Size Markers

Ambion offers several different ranges of RNA size markers that can be obtained unlabeled for staining with EtBr or biotinylated for subsequent secondary detection (Cat. No. 7140 - unlabeled, #7175 - biotinylated) contains 5 transcripts evenly spaced between 100 -500 nt, which are ideal for ribonuclease protection assays. RNA Century Markers can also be obtained as DNA templates (Cat. No. 7780 and 7782) for the synthesis of radiolabeled RNA markers in an in vitro transcription. **Marker Set** (Cat. No.7150 - unlabeled, #7170 - biotinylated) contains 10 transcripts ranging from 0.5-9.0 kb for use with Northern analysis.

RNA transcripts and double-stranded DNA markers (e.g. pUC 19/Hpa II, Cat. No. **7760** and **7770**) can also be end-labeled with polynucleotide kinase (5' end-labeling reaction) and denatured, for use as labeled size markers.

Other guides to RNA size and migration position are the xylene cyanol and bromophenol blue dyes present in most loading buffers, and rRNA species present in Northern analysis. The migration position of the dyes included in loading buffers is affected both by gel percentage and composition (denaturing vs. nondenaturing) and total RNA samples. Both the 18S and 28S species are strongly visible in Northern gels stained with EtBr or UV-shadowed. The table above gives their sizes in

## References

- Chomczynski, P. (1992) Solubilization in formamide protects RNA from degradation. Nuc. Acids Res. 20:3791-3792.
- Wallace, D.M. (1987) Precipitation of Nucleic Acids. Methods of Enzymology 152:41-46.

## Ordering Information

製品番号 (カタログ番号)	AM1964
製品名	RNaseAlert™ Lab Test Kit
大小	25 reactions
価格 (JPY)	40,100 (お問い合わせください)
量	-

製品番号 (カタログ番号)	AM1966
製品名	RNaseAlert™ QC System
大小	480 assays
価格 (JPY)	171,900 (お問い合わせください)
量	-

製品番号（カタログ番号）**AM7000**  
製品名**THE RNA Storage Solution**  
大小10 x 1 mL  
価格 (JPY)10,000  
(  
お問い合わせください。)  
量-

製品番号（カタログ番号）**AM7001**  
製品名**THE RNA Storage Solution**  
大小50 mL  
価格 (JPY)10,000  
(  
お問い合わせください。)  
量-

製品番号（カタログ番号）**AM7145**  
製品名**Century™-Plus RNA Markers**  
大小1 tube  
価格 (JPY)32,900  
(  
お問い合わせください。)  
量-

製品番号（カタログ番号）**AM7150**  
製品名**Millennium™ RNA Markers**  
大小50 µL  
価格 (JPY)39,700  
(  
お問い合わせください。)  
量-

製品番号（カタログ番号）	AM7760
製品名	pUC19 DNA (Sau3A I digested)
大小	50 µg
価格 (JPY)	41,700 ( お問い合わせください。)
量	-

製品番号（カタログ番号）	AM7780
製品名	RNA Century™ Marker Templates
大小	5 µg
価格 (JPY)	41,700 ( お問い合わせください。)
量	-

製品番号（カタログ番号）	AM7782
製品名	RNA Century™-Plus Marker Templates
大小	5 µg
価格 (JPY)	41,700 ( お問い合わせください。)
量	-

製品番号（カタログ番号）	AM9480
製品名	LiCl Precipitation Solution (7.5 M)
大小	100 mL
価格 (JPY)	19,500 ( お問い合わせください。)
量	-

製品番号（カタログ番号）	AM9780
製品名	RNaseZap™ RNase Decontamination Solution
大小	250 mL
価格 (JPY)	14,500 ( お問い合わせください )
量	-

製品番号（カタログ番号）	AM9782
製品名	RNaseZap™ RNase Decontamination Solution
大小	6 x 250 mL
価格 (JPY)	64,500 ( お問い合わせください )
量	-

製品番号（カタログ番号）	AM9784
製品名	RNaseZap™ RNase Decontamination Solution
大小	4 L
価格 (JPY)	127,900 ( お問い合わせください )
量	-

製品番号（カタログ番号）	AM9786
製品名	RNaseZap™ RNase Decontamination Wipes
大小	100 sheets
価格 (JPY)	14,500 ( お問い合わせください )
量	-

製品番号（カタログ番号）	AM9788
製品名	RNaseZap™ RNase Decontamination Wipes Refill
大小	300 sheets
価格 (JPY)	38,700 ( お問い合わせください。)
量	-

製品番号（カタログ番号）	AM9860
製品名	TE, pH 7.0, RNase-free
大小	10 x 1 mL
価格 (JPY)	10,000 ( お問い合わせください。)
量	-

製品番号（カタログ番号）	AM9861
製品名	TE, pH 7.0, RNase-free
大小	50 mL
価格 (JPY)	10,000 ( お問い合わせください。)
量	-

製品番号（カタログ番号）	AM9912
製品名	EDTA (0.1 mM), pH 8.0, RNase-free
大小	50 mL
価格 (JPY)	9,800 ( お問い合わせください。)
量	-

For Research Use Only. Not for use in diagnostic procedures.

Japan | 日本



Select Language:

Japan - 日本語 | Japanese

▼

Go

[ウェブサイト利用条件](#) | [個人情報保護方針](#) | [標準販売条件](#)